

V3

Featuring Patented AIR°gel Technology

The V3 combines AIR°gel technology with acute care performance for a revolutionary home care experience. The V3 fits quickly and easily while enhancing patient comfort.

AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

AIR°gel
Forehead Cushion

Forehead Adjustment Screw

Quick-Release Slot

Lightweight, breathable headgear slides off easily with one hand.


360° Rotating Elbow

For optimal tube positioning and easier mobility

Oxygen Port

Allows easy delivery of supplemental oxygen.

Size Indicator
Available in four sizes



Made in the USA

What makes our V3 Mask different?

- ✓ Adjustment Screw -provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the bridge of the nose.
- ✓ AIR°gel Forehead Pad - stabilizes the mask comfortably on the patient's face during a night's sleep minimizing pressure points.
- ✓ Available in four sizes with simple and breathable quick-release headgear.

Advanced Cushion Technology(A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



V3 Masks with EZfit Headgear For Comfortable & Secure Treatment at Night



V3 Product Specifications:

Product Number	Mask Style	Leak Port	Entrianment Valve	Size	Mask Weight	Mask Dimensions
50863	Vented	✓	✓	Small	109.2g	3.60"D x 4.10"W x 7.25"H
50864	Vented	✓	✓	Medium	114.3g	3.60"D x 4.30"W x 7.40"H
50865	Vented	✓	✓	Large	119.4g	3.60"D x 4.38"W x 7.60"H
50866	Vented	✓	✓	Extra Large	120.1g	3.60"D x 4.42"W x 7.80"H