

Stimulate skin rejuvenation and slow down the aging process



Made in Canada




Consult the Official Warranty Statement issue by the manufacturer to fully understand the conditions.

Benefits

- ✓ Achieve better blood circulation
- ✓ Eliminate toxins efficiently
- ✓ Encourage neuroplasticity (regeneration of brain cells)
- ✓ Increase various functions of metabolism anti-aging effects
- ✓ Repair gray matter (which contains nerve cell bodies)
- ✓ Repair white matter (which are the filaments that connect cell bodies and carry messages between nerve cells)
- ✓ Reach all oxygen-deficient areas of the brain to reactivate neurones
- ✓ Improve motor skills, muscle tone, motricity, spasticity
- ✓ Increase fine motor skills
- ✓ Reduce muscle cramps Increase strength
- ✓ Increase oxygen saturation to the brain
- ✓ Reduce cerebral edema and neuroinflammation
- ✓ Grow new blood vessels (Angiogenesis)
- ✓ Increase stem cell production (of brain cells)

Oxynova 9 includes



Air Compressor
Lab-grade high flow



Filters
5- micron filter
(optional .01 micron)



Anti- Roll Bolsters
Two full length
extra-rigid stabilizers



Carry Bag
Double- layered nylon




Fabric (Bladder)
PVDF Teflon- coated; assembly
method: double overlap 3.175 cm
thermal fusion fabric welding and
radio frequency welding




Cover Assembly
Adds a layer of structure to
the bladder (unique to our
chambers), easy to clean and
water repellent



Hose
Laboratory-grade 9 ft (2.5 m) hose,
with reinforced nylon sheath



**Interior & Exterior
Air Pressure Gauge**
100% stainless steel




Safety Feature
2 custom pressure-relief valves;
pressure adjustment handle

Technical specs

Working Pressure	1.3 ATA (0.3 bar / 40 kPa / 4.3 psi)
Windows	10 double layared 14 cm (5 ½ inches) windows
Air Volume Delivery	160 Litres/min
Volume	0.562 m3 (± 5 %)
Diameter	Ø 92 cm (36 inches)
Lenght	275 cm (108 inches)
Chamber's Weight	28 kg / 51.7 lbs
Compressor's Weight	18.1 kg / 40 lbs
Sound level	Less than 60 db (A)
Power	800 W, Single Phase
Pass-through fitting	10 custom-interchangeable positions
Single mask adaptor	Oxygen adaptor mask


Optional equipment



(optional)
"ICE PACK" air collar



(optional)
Oxygen nasal cannula



(optional)
Oxygen mask



(optional)
Oxygen concentrator
(5 L/m)



(optional)
Oxygen concentrator
(10 L/m)



Plot No.20, KVT Garden, Phase-III, Madavaram, Chennai - 60. TN
Branches: Chennai | Bangalore | Hyderabad | Cochin



Contactus@ababilhealthcare.com
www.ababilhealthcare.com



+91 44 4858 2255
+91 98840 43905