

Handheld Patient Monitors M850 (SPO2+NIBP)



Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: ECG, RR(Standard), SpO2, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Selectable 3-lead and 5-lead for ECG measurement.
- ✓ Selectable layout modes including font and waveforms.
- ✓ Lithium battery support more than 28 hours continuous working.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M860 (SPO2+ECG)

Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: NIBP, RR(Standard), SpO2, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Specifically designed for a reliable and accurate blood pressure measurement.
- ✓ Lithium battery support more than 20 hours continuous working.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M880 (SPO2+ETCO2)



Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: CO₂, RR(Standard), SpO₂, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Minimal warm-up time for EtCO₂ and RR measurements.
- ✓ CO₂ water filter ensuring extended monitoring in a humid environment.
- ✓ Selectable layout modes including font and waveforms.
- ✓ More than 18 hours lithium battery.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M800

Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter:
SpO₂, PR(Standard), Nellcor SpO₂(Optional),
ECG, RR(Optional).
- ✓ 2.4" color TFT display with auto rotation function.
- ✓ Selectable layout modes including font and waveforms.
- ✓ Dual work modes: continuous monitoring & spot check.
- ✓ 3x AA batteries or rechargeable battery with charger stand(optional).



Ajustable alarms



Power saving



Easy to carry



300 Groups of Memory

Oximeter M70



Care about your SpO₂, care about your health.

- ✓ Dual color OLED display.
- ✓ Parameter: SpO₂, PR, waveform and pulse bar.
- ✓ 2 directions display.
- ✓ Resistance to surrounding light and man-made light.
- ✓ Colors optional: Yellow, Green, White.
- ✓ Auto shut off without signal after 8s.

M70A



M70B



M70D



Low power consumption



5s for measuring



Easy to carry



Medical level

Oximeter M70C

Care about your SpO₂, care about your health.

- ✓ OLED screen with 6 display modes.
- ✓ Colors optional: blue, Yellow.
- ✓ 5 levels adjustable brightness.
- ✓ Resistance to surrounding light and man-made light.
- ✓ Auto shut off without signal after 8s.
- ✓ With Bluetooth and professional APP optional.



Bluetooth4.0 (Optional)



Low power consumption



5s for measuring



Easy to carry



Medical level


Fetal Doppler WF100



CE

Hear your baby from the beginning

- ✓ Compact and portable design.
- ✓ Large color LCD display.
- ✓ High sensitive probe.
- ✓ Built-in speaker and interface of earphone.
- ✓ Three working modes: Real-time, Average and Manual counting.
- ✓ 3 x AA batteries or rechargeable Lithium battery(optional).



BP Tracker WBP302

Simpler and faster

- ✓ Blood pressure, HR, Blood range display.
- ✓ Large & eye protection backlight LCD screen.
- ✓ Intelligent pressurized.
- ✓ Clinically accurate & fast reading.
- ✓ Indication for wrong operation.
- ✓ Special design to place cuff.
- ✓ 100 groups of memory.
- ✓ Various cuff options.



CE

