

# Aura

featuring AIR<sup>o</sup>gel with Advanced Cushion Technology


Enjoy comfort and a great seal with the world's first nasal mask featuring both Custom Fit Technology<sup>TM</sup> and an Active Headgear Connector<sup>TM</sup>. With three sizes and a Flexible Spacebar<sup>TM</sup>, Aura will give you the rest you deserve.

## Enter a New Aura of Comfort and Seal

### Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

**Flex, fit, forget.**



### AIR<sup>o</sup>gel

Maximum comfort, minimum weight. Our patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation.

### Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

### 360° Rotating Elbow

For optimal tube positioning and easier mobility

### Flexible Spacebar

Unique bendable design provides superior stability, holding a seal without contacting the skin.

### Size Indicator

available in three sizes.




### Advanced Cushion Technology<sup>TM</sup> (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

### Active Headgear Connector<sup>TM</sup>

Moves as the patient moves for a reliable seal.



Made in the USA

# What makes our Aura Mask different?

- ✓ Custom Fit Technology is unique to Sleepnet. Our soft-shell mask has a shapeable wire allowing for optimum flexibility, fit, and comfort.
- ✓ Our Active Headgear Connector moves when you move to help maintain a seal throughout the night.
- ✓ Flexible Spacebar eliminates any contact with the forehead, yet retains stability.
- ✓ Available in three sizes with simple and breathable quick-release EZ Fit headgear.

## Advanced Cushion Technology(A.C.T.)

Our patented AIR<sup>o</sup>gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



## Custom Fit Technology



**BEND IT.**



**MOLD IT.**



**SHAPE IT.**

## Aura Mask with Headgear Specifications:

Product Number	Case Quantity	Style	Leak Port	Size	Leak Port	Mask Dimensions
50176	50	Vented	✓	S	61.5g	5.50"D x 2.85"W x 3.80"H
50177	50	Vented	✓	M	62.5g	5.50"D x 2.85"W x 3.80"H
50178	50	Vented	✓	L	71.5g	5.50"D x 2.85"W x 3.80"H
55051	10	EZFit Headgear	N/A	N/A	N/A	N/A