

Handheld Patient Monitors M850 (ECG, SPO2 and RESP)



Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: ECG, RR(Standard), SpO2, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Selectable 3-lead and 5-lead for ECG measurement.
- ✓ Selectable layout modes including font and waveforms.
- ✓ Lithium battery support more than 28 hours continuous working.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M860 (NIBP, SPO2 and RESP)

Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: NIBP, RR(Standard), SpO2, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Specifically designed for a reliable and accurate blood pressure measurement.
- ✓ Lithium battery support more than 20 hours continuous working.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M880 (ETCO2, SPO2 and RESP)



Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter: CO2, RR(Standard), SpO2, PR(Optional).
- ✓ 4.3" color TFT display with auto rotation function.
- ✓ Minimal warm-up time for EtCO2 and RR measurements.
- ✓ CO2 water filter ensuring extended monitoring in a humid environment.
- ✓ Selectable layout modes including font and waveforms.
- ✓ More than 18 hours lithium battery.



Dual work modes



Power saving



Ajustable alarms



Easy to carry



500 Groups of Memory

Handheld Patient Monitors M800

Easier to carry, easier to operate

- ✓ Suitable for adult, pediatric and neonate.
- ✓ Apply to hospital, clinic and home.
- ✓ Monitoring parameter:
SpO2, PR(Standard), Nellcor SpO2(Optional),
ECG, RR(Optional).
- ✓ 2.4" color TFT display with auto rotation function.
- ✓ Selectable layout modes including font and waveforms.
- ✓ Dual work modes: continuous monitoring & spot check.
- ✓ 3x AA batteries or rechargeable battery with
charger stand(optional).



Ajustable alarms



Power saving



Easy to carry



300 Groups of Memory

Oximeter M70



Care about your SpO₂, care about your health.

- ✓ Dual color OLED display.
- ✓ Parameter: SpO₂, PR, waveform and pulse bar.
- ✓ 2 directions display.
- ✓ Resistance to surrounding light and man-made light.
- ✓ Colors optional: Yellow, Green, White.
- ✓ Auto shut off without signal after 8s.

M70A



M70B



M70D



Low power consumption



5s for measuring



Easy to carry



Medical level

Oximeter M70C

Care about your SpO₂, care about your health.

- ✓ OLED screen with 6 display modes.
- ✓ Colors optional: blue, Yellow.
- ✓ 5 levels adjustable brightness.
- ✓ Resistance to surrounding light and man-made light.
- ✓ Auto shut off without signal after 8s.
- ✓ With Bluetooth and professional APP optional.



Bluetooth4.0 (Optional)



Low power consumption



5s for measuring



Easy to carry



Medical level


Fetal Doppler WF100



Hear your baby from the beginning

- ✓ Compact and portable design.
- ✓ Large color LCD display.
- ✓ High sensitive probe.
- ✓ Built-in speaker and interface of earphone.
- ✓ Three working modes: Real-time, Average and Manual counting.
- ✓ 3 x AA batteries or rechargeable Lithium battery(optional).

CE



BP Tracker WBP302

Simpler and faster

- ✓ Blood pressure, HR, Blood range display.
- ✓ Large & eye protection backlight LCD screen.
- ✓ Intelligent pressurized.
- ✓ Clinically accurate & fast reading.
- ✓ Indication for wrong operation.
- ✓ Special design to place cuff.
- ✓ 100 groups of memory.
- ✓ Various cuff options.



CE

