

HP-TCI Anesthetist's Best Partner



Advanced TCI function is offering a great anesthesia solution

- ✓ Clear visibility on important information with medical grade large screen
- ✓ Easy operation with 3.0 inch touchscreen
- ✓ Accurate concentration is controlled by reliable algorithm
- ✓ Extra bolus delivery is available during the TCI
- ✓ Max flow rate: 2000mL/h


Overview of Possible Selections

Drug	Model	Plasma	Effect side
Propofol	Marsh	✓	X
	Schnider	✓	✓
	Kataria (Optional)	✓	X
	Paedfusor (Optional)	✓	X
Remifentanil	Minto	✓	✓
Sufentanil	Gepts (Optional)	✓	✓
Alfentanil	Maitre (Optional)	✓	✓
Fentanyl	Shafer (Optional)	✓	✓




ABABIL
PEOPLE | PASSION | PRODUCTS


As a member of HP family,
all the advantages of HP family is inherited.



Syringe position detection system.



Full range syringe compatibility




Auto/Manual installation system




It is not only a pump, but also a Station:

Easy to increase to multi-Channels via
HP-80 station for the use of TIVA.

Comprehensive internet connection solution.

 contactus@ababilhealthcare.com
www.ababilhealthcare.com

 20, KVT Garden Phase-III, VS Mani Nagar, Madhavaram, Chennai - 60. TN
Branches: Chennai | Bangalore | Hyderabad | Cochin

 +91 44 2553 2144
+91 98840 43905

DRIVEN BY PEOPLE, FUELLED BY PASSION, INSPIRED BY PRODUCTS