

iQ Blue

featuring AIR^ogel with Advanced Cushion Technology

Introducing the next generation in Airgel. Our iQ[®] Blue Nasal Vented and Non-Vented masks feature our patented AIR^ogel with Advanced Cushion Technology and flexible shell with Custom Fit Technology[™]. Together, they provide an optimal fit with maximum comfort. Unlike masks with rigid shells, iQ Blue adjusts to truly fit the contours of your patients' face.

The Smart Choice for a Custom Fit



Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

Flex, fit, forget.

Clearest Field of Vision Possible

No forehead pad required.

AIR^ogel

Maximum comfort, minimum weight. Our patented AIR^ogel minimizes pressure points and skin irritation.

Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.


Advanced Cushion Technology[™] (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility

Vented and Non-Vented Versions Available



Made in the USA

What makes our iQ Blue Mask different?

- ✓ Custom Fit Technology - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- ✓ Our AIR°gel with Advanced Cushion Technology creates an ergonomic gel cushion that reduces pressure points on the skin.
- ✓ Clearest field of vision possible - users can even wear glasses with this mask.
- ✓ No forehead pad, yet retains outstanding stability.

Advanced Cushion Technology(A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Custom Fit Technology



BEND IT.



MOLD IT.



SHAPE IT.

iQ Blue Mask Specifications:

Product Number	Mask Style	Leak Port	Headgear Style	Mask Size	Tube Length	Mask Weight	Mask Dimensions
50647	Vented	✓	3-Point	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50655	Vented	✓	StableFit	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50250	Vented	✓	3-Point	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50255	Vented	✓	StableFit	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50660	Non-Vented	✗	3-Point	OS	6"	63.8 g	3.08"D x 3.13"W x 2.50"H