

Innova Hospital


featuring AIR^ogel with Advanced Cushion Technology

The Innova[®] Hospital disposable mask features our patented AIR^ogel with Advanced Cushion Technology and Active Headgear Connector[™]. Together, they provide an optimal fit, maximum comfort, and freedom to move without breaking the seal.

Move Freely with Our Disposable Hospital Nasal Mask

AIR^ogel

Maximum comfort, minimum weight. Our patented AIR^ogel minimizes pressure points and skin irritation.



Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

360° Rotating Elbow

For optimal tube positioning and easier mobility

Touchless Spacebar[™]

Unique design eliminates any contact with the patient's forehead.

Size Indicator

available in two sizes.

S/M M/L

Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

Active Headgear Connector

Moves as you move for a reliable seal.

Vented and Non-Vented Versions Available

For CPAP, Bi-level, and ventilator therapy



Made in the USA

What makes our Innova Mask different?

- ✓ Vented and nonvented versions are available for easy hook up to CPAP, Bi-level, or ventilaor therapy..
- ✓ Easy-to-use with a Quick-Release Headgear Connector.
- ✓ Touchless Spacebar™ eliminates any contact with the forehead, yet retains stability.
- ✓ Available in two sizes with 4-point simple and breathable headgear.

Advanced Cushion Technology (A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Active Headgear Connector

Fifteen degrees of rotation provide relief from headgear tension, allowing the mask to move with the patient for a reliable seal.

Innova Hospital Mask with Headgear Specifications:

Product Number	Style	Leak Port	Size	Mask Weight	Mask Dimensions
50180	Vented	✓	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50181	Vented	✓	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H
50594	Non-Vented	X	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50595	Non-Vented	X	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H