

Technical specifications

Performance

Default settings 25 cmH₂O (target pressure), 5 cmH₂O above set pressure for 5 min (time-limited hold)

Physical dimensions

Dimensions (L x W x H) 15.7 cm (6.2 in) x 5.5 cm (2.2 in) x 3.6 cm (1.4 in)
 Weight (incl. batteries) 260 g (9.2 oz)
 Display 60 x 35 mm
 Ingress Protection (IP) IP 34

Electrical characteristics

AC power input 100 to 240 VAC / 50 to 60 Hz \pm 10%; 1.25 VA typical, 3.25 VA maximum
 USB input power 5 VDC \pm 5%; 1.0 W typical, 2.5 W maximum
 Car adapter 12 V to 24 V \pm 5%; 1.25 W typical, 3.25 W maximum
 Battery requirements AA (IEC-HR6) NiMH rechargeable 1.2 V, \geq 1900 mA

Connectors

Electrical Mini USB
 Pneumatic Male luer slip

Environment

Temperature -15°C to 50°C / 5°F to 122°F (operating), -15°C to 70°C / 5°C to 158°F (storage)
 Relative humidity 5% to 95% noncondensing (operating) / 5% to 95% noncondensing (storage)
 Noise < 30 dB(A)
 Altitude -650 to 7,620 m / -2,132 to 25,000 ft above sea level
 Above 4,000 m / 13,123 ft, only battery power is supported

