

# Mojo 2

## Featuring Patented AIR<sup>o</sup>gel Technology

The Mojo 2 Vented and Non-Vented masks feature a patented AIR<sup>o</sup>gel cushion, moldable shell with Custom Fit Technology™ and a Forehead Adjustment Screw. Together, they provide a personalized fit with maximum comfort.

### AIR<sup>o</sup>gel

Maximum comfort, minimum weight. Our patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation.



### Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

**Bend It. Mold It. Shape It.**

### Oxygen Port

Allows easy delivery of supplemental oxygen.

### 360° Rotating Elbow

For optimal tube positioning and easier mobility

### Forehead

#### Adjustment Screw

Removes pressure from the bridge of the nose for personalized fit.

### AIR<sup>o</sup>gel

Forehead Cushion


Vented, Non-Vented, and Non-Vented AAV Versions Available

### Quick-Release Magnetic Headgear Clips

Lightweight, breathable headgear disconnects from mask easily with one hand.

### Size Indicator

Available in four sizes



Made in the USA

# What makes our Mojo 2 Mask different?

- ✓ Custom Fit Technology - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- ✓ Adjustment Screw - provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the forehead.
- ✓ Available in four sizes with magnetic connectors and simple, breathable quickrelease headgear.

## About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



## Custom Fit Technology



**FLEX**



**FIT**



**FORGET**

## Mojo 2 Masks with Headgear Specifications:

Product Number	Style	Size	Leak Port	Entrianment Valve	Mask Weight	Mask Dimensions
50817	Vented	Small	✓	✓	133.8g	3.60"D x 4.30"W x 7.25"H
50818	Vented	Medium	✓	✓	141g	3.60"D x 4.40"W x 7.40"H
50819	Vented	Large	✓	✓	144.5g	3.60"D x 4.50"W x 7.60"H
50820	Vented	Extra Large	✓	✓	146.5 g	3.60"D x 4.55"W x 7.80"H
50912	Non-Vented	Small	X	X	133.2g	3.70"D x 4.30"W x 7.25"H
50913	Non-Vented	Medium	X	X	140.4g	3.70"D x 4.40"W x 7.40"H
50914	Non-Vented	Large	X	X	143.9g	3.70"D x 4.50"W x 7.60"H
50915	Non-Vented	Extra Large	X	X	145.9 g	3.70"D x 4.55"W x 7.80"H
50941	AAV Non-Vented	Small	X	✓	133.7g	3.65"D x 4.30"W x 7.25"H
50942	AAV Non-Vented	Medium	X	✓	140.9g	3.65"D x 4.40"W x 7.40"H
50943	AAV Non-Vented	Large	X	✓	144.4g	3.65"D x 4.50"W x 7.60"H
50944	AAV Non-Vented	Extra Large	X	✓	146.4g	3.65"D x 4.55"W x 7.80"H