

HBOT is a foundational treatment to stimulate the brain to repair itself.



Made in Canada



LIFETIME warranty



Consult the Official Warranty Statement issue by the manufacturer to fully understand the conditions.


Benefits

- ✓ Achieve better blood circulation
- ✓ Eliminate toxins efficiently
- ✓ Encourage neuroplasticity (regeneration of brain cells)
- ✓ Increase various functions of metabolism anti-aging effects
- ✓ Increase control of the body and mind
- ✓ Decrease edema effects by vasoconstriction
- ✓ Decrease leucocyte chemotaxis and adhesion
- ✓ Attenuate ischemic-reperfusion damage
- ✓ Improve motor skills in children, muscle tone, motricity, spasticity
- ✓ Increase fine motor skills
- ✓ Reduce muscle cramps
- ✓ Increase strength
- ✓ Suppress the formation of inflammatory mediators
- ✓ Inhibit the autoimmune syndrome and the immune reaction in antigens
- ✓ Improve chronic skin damage healing by inducing angiogenesis
- ✓ Stimulate O2 dependent collagen matrix formation, which is an essential phase in wound healing




ABABIL
PEOPLE | PASSION | PRODUCTS

Oxynova 7 includes



Air Compressor
Lab-grade high flow



Filters
5- micron filter
(optional .01 micron)



Anti- Roll Bolsters
Two full length
extra-rigid stabilizers



Carry Bag
Double- layered nylon




Fabric (Bladder)
PVDF Teflon- coated; assembly
method: double overlap 3.175 cm
thermal fusion fabric welding and
radio frequency welding




Cover Assembly
Adds a layer of structure to
the bladder (unique to our
chambers), easy to clean and
water repellent



Hose
Laboratory-grade 9 ft (2.5 m) hose,
with reinforced nylon sheath



**Interior & Exterior
Air Pressure Gauge**
100% stainless steel




Safety Feature
2 custom pressure-relief valves;
pressure adjustment handle

Technical specs

Working Pressure	1.4 ATA (0.4 bar / 40 kPa / 5.7 psi)
Windows	8 double layared 14 cm (5 ½ inches) windows
Air Volume Delivery	160 Litres/min
Volume	0.562 m3 (± 5 %)
Diameter	Ø 72 cm (28 inches)
Lenght	229 cm (90 inches)
Chamber's Weight	14.5 kg / 31.9 lbs
Compressor's Weight	18.1 kg / 40 lbs
Sound level	Less than 60 db (A)
Power	800 W, Single Phase
Pass-through fitting	8 custom-interchangeable positions
Single mask adaptor	Oxygen adaptor mask

Optional equipment



(optional)
"ICE PACK" air collar



(optional)
Oxygen nasal cannula



(optional)
Oxygen mask



(optional)
Oxygen concentrator
(5 L/m)



(optional)
Oxygen concentrator
(10 L/m)