

# V3

## Featuring Patented AIR°gel Technology

The V3 combines AIR°gel technology with acute care performance for a revolutionary home care experience. The V3 fits quickly and easily while enhancing patient comfort.

### AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

AIR°gel  
Forehead Cushion

Forehead Adjustment Screw

### Quick-Release Slot

Lightweight, breathable headgear slides off easily with one hand.

360° Rotating Elbow  
For optimal tube positioning and easier mobility

### Oxygen Port

Allows easy delivery of supplemental oxygen.

Size Indicator  
Available in four sizes



Made in the USA

## What makes our V3 Mask different?

- ✓ Adjustment Screw -provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the bridge of the nose.
- ✓ AIR°gel Forehead Pad - stabilizes the mask comfortably on the patient's face during a night's sleep minimizing pressure points.
- ✓ Available in four sizes with simple and breathable quick-release headgear.

### Advanced Cushion Technology(A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



## V3 Masks with EZfit Headgear For Comfortable & Secure Treatment at Night



## V3 Product Specifications:

Product Number	Mask Style	Leak Port	Entrianment Valve	Size	Mask Weight	Mask Dimensions
50863	Vented	✓	✓	Small	109.2g	3.60"D x 4.10"W x 7.25"H
50864	Vented	✓	✓	Medium	114.3g	3.60"D x 4.30"W x 7.40"H
50865	Vented	✓	✓	Large	119.4g	3.60"D x 4.38"W x 7.60"H
50866	Vented	✓	✓	Extra Large	120.1g	3.60"D x 4.42"W x 7.80"H