

SunTech CT40

With more than 30 years dedicated to developing our world-class blood pressure technology, we are confident that we don't miss a beat - so you won't either. It's our legacy. It's SunTech Medical. It's The Difference In Clinical Grade.

Reliability:

Manual SphygMode™ for real-time verification of BP readings without removing the cuff

Accuracy:

Customizable BP averaging protocol to ensure accurate, consistent readings

Connectivity:

Transfer of data using IHE protocols with network connection through built-in ethernet or optional Wi-Fi

Flexibility:

Modular design for addition of optional temperature and SpO2 modules

Simplicity:

Intuitive user interface with LCD display

Modular Vital Signs with BP, Temperature and SpO₂



Standard Features

- ✓ Clinical-grade Advantage™ BP technology
- ✓ Customizable BP averaging (SPRINT recommended)
- ✓ Manual SphygMode™ for real-time BP verification
- ✓ HL7 messaging for data transfer
- ✓ Built-in memory - up to 99 charting sessions
- ✓ LCD display with simple design for ease of use
- ✓ LAN/Ethernet & 2 USB ports

Options & Accessories

- | | |
|------------------------------------|-------------------|
| ✓ Touchless IR Thermometry* | ✓ Barcode Scanner |
| ✓ Covidien Filac® 3000 Thermometry | ✓ Mobile stand |
| ✓ Masimo SET® SpO2 (preferred) | ✓ Printer |
| ✓ ChipOx Nellcor® -compatible SpO2 | ✓ Wi-Fi |

SUNTECH