



Vein Navi

Simple Visibility

NAVI-30 Vein Illuminator provides a real time accurate image of the vasculature on the surface of the skin, which helps clinicians to verify the status of the vein and reduces the risk of venipuncture.

With the help of NAVI-30 Vein Illuminator, hospital may gain the benefit of reduce costs for healthcare, improve infection control and improve patient safety and satisfaction.

Portable

Hand hold type Vein Illuminator, easy to carry, non-contact to patient.

Colorful

Great image with 5 different colors optional, suitable for different types of skin.





Specifications

Resolution	854 x 480 pixel
Image Light	10-25 Lumens 4 level adjustable
Color	5 colors
Image Size	3 Sizes
Charge during operation	Yes
Size	224 (W) x 68 (H) x 64 (D) mm
Weight	500 g

Application



Pediatrics

Small veins "Baby fat" Anxiety of the children and their families.



Oncology

Fewer options for successful access because of repeated infusions and procedures and medication.



Vascular Surgery

Sclerotherapy for treating spider telangiectasia Vein stripping treating.



Plastic Surgery

Sclerotherapy for treating spider telangiectasia Cosmetic injections on the face