

S Series

CPAP-S1 & BIPAP-S2



Precious and Elegance



5.0 inch enlarged screen design enable a better parameter review.



Integrated Humidification tank design to improve the efficiency of heating and humidifying



360 ° adjustable air outlet, suitable for various sleeping positions



Intuitive user interface

Smart and Simplify

Auto-fix algorithm

Inheriting the intelligent medical algorithm, the rhythm information model detection method established by tracking and learning the characteristics of the patient's change, learns the patient's change and continuously and automatically corrects it, accurately predicts the patient's trigger point, reduces the patient's confrontation, and increases the comfort of the patient.



P-Breath Home Breathing Cloud

Treatment data is available at your fingertips, meeting remote analysis and data intelligent management requirements.



Wireless Spo2 module Interconnection

- ✓ Optional wireless Spo2 module for more comprehensive treatment



Low Noise

- ✓ Operating under low noise for a more comfortable interactive experience



High Performance turbine

- ✓ Faster response and more stable delivery airflow



- ✓ Automatic Ramp
- ✓ Intelligent heating
- ✓ Smart start and stop
- ✓ Automatic screen off
- ✓ Automatic leak Compensation



Products Parameters :

Specification	S Series	
Product Name	S1	S2
Mode	CPAP Auto-CPAP	CPAP S/T S T VCPS
Dimensions	270mm*135mm*100mm(L*W*H)	
Weight	1.5kg	
Pressure	4-20cmH2O	4-30cmH2O
Efficient humidifier	Yes	Yes
Smart Star/Stop	Yes	Yes
EPR level	OFF, 1-3	

Specification	S Series	
Spo2 module	Yes	Yes
Leak Compensation	Yes	Yes
Heating tube	Optional	Optional
LCD screen brightness	1-3 Brightness setting, Auto screen off	
Build in wireless technology	Yes	Yes
Bluetooth, can be used to pair the pulse oximeters.	Optional	Yes